

## Clinical Upskill Workshop for Biomed teams

Patient Monitoring | Classroom - Virtual

### Synopsis

Upskill your technicians with relevant knowledge and skills to enhance their understanding of Clinical applications within your PHILIPS patient monitoring solution. Learn more about advanced features, additional troubleshooting techniques and clinical decision support tools.

### Class capacity

The course supports a limited number of participants. Attendance is on a first come-first serve basis.

### Target audience

- Biomedical Engineers
- Technicians

### Course details

- Duration: 4 hours
- Virtual: A link will be shared with you closer to the date
- Cost: AUD \$495.00 pp (excl GST)

### Speaker



### Karyn Harris

Karyn is a Philips Clinical Consultant and has worked within the Patient Monitoring division for the past 17 years.

A Registered Nurse with a long and lively history in Critical Care nursing in Australia and overseas, Karyn has an extensive knowledge base as well as real life scenarios to draw on.

Karyn provides an easy-going learning environment and is passionate about providing skills and information to Biomed and technicians to enable them to communicate more easily in complex Critical care environments and to get more out of their Intellivue Monitoring systems.

### Agenda overview

- Advanced measurements
- Clinical Application
- Admission/Discharge/Transfer
- Audit log
- Basic A & P
- Change alarms
- Modify/Manage configs
- Intro to screen design
- Profiles
- ECG arrhythmia
- Alarm Advisor

### What will you get?

- Certificate of attendance
- 1 year access to Incentre